

Bottle rack, plastic-coated, 20 compartments

Produto n.º: 2551801
EUR Preço: Contacte-nos
Enviado em 2 semanas

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Especificações

Volume / Tamanho da embalagem: 1